

Universal Forefoot Compressive Cryotherapy Device

DESIGNED TO

- Providing dynamic support
- Fighting against pain
- Limiting bleeding and hematoma
- Managing post operative edema
- Easy to use at home



Excell' ICE



Excell'ICE cold kits combine cryotherapy with modulated static compression using an inflation pump.

Easy to apply, its technology consists of alternating compressive cold during recovery and/or healing sessions, commonly recommended in the days following an injury or post-operation.

It is well known that external static compression with ice increases the extent of cooling at the surface of the skin.

A unique new concept where the compression chamber is directly integrated into the cold pack, making it much easier to use. This new technology makes it easier to carry out 20-30 minute cryotherapy-compressive recovery sessions. These portable devices are very practical to use at home or when travelling.

For athletes, this will ease musculoskeletal pain immediately after their matches or competitions and enable them to recover optimally, thereby reducing the time they spend out of competition due to pain, injury, fatigue, etc.



Recommended use:

20 min. by session every 4 hours
First session in Operating Room after dressing

Made in Europe

INDICATIONS

- Post-operative care
- Analgesic treatment
- Sport trauma (sprain, tear, etc..)
- After rehabilitation session

CONTRAINDICATIONS

- Allergy to cold
- Insensitivity to cold
- Raynaud's syndrome
- Cryoglobulinemia
- Arteritis
- Diabetes
- Polyarthritis...

Excell' ICE

Today's most user-friendly system: "Cold Pack" and "Compression Chamber" combined in a single device.

- Compression multiplies the effects of cold
- Cold packs remain flexible to adapt to body shape even at -4°F (-20°C)
- No frost to avoid skin burns
- Hypoallergenic Lycra with tightened mesh (included)
- 360° coverage around the Wrist in a single cold pack
- Freezer bag with zip slider for easy opening/closing (included)

Forefoot Cryo kit

Introducing the Forefoot Cryo Kit, a cooling solution that is incredibly easy to apply with targeted cooling effect on the Forefoot (not the ankle nor the tip of the toes). With its universal one-size-fits-all design, it will fit any foot. This kit includes one single cold pack surrounding the forefoot and easy Velcro straps, providing comprehensive coverage for the forefoot. The innovative design combines cold therapy with compression in a single, convenient piece. Achieve optimal compression using the included manual pump.



1 Place the cold pack on a flat surface. The rounded part should be level with the heel.



2 Fold the left part over the foot.



3 Pull the right side over the left and fasten with velcro.



4 Position and fasten the heel strap.



5 Set the valve to 0 and screw on the pump. Inflate (do not exceed 10 pressures). Set valve to F and unscrew pump.

Universal Forefoot Compressive Cryotherapy Device Ref. ORT40PACK-P



Cold pack



Manual pump



Heel strap



Jersey



Freezer bag



Suitcase