

KHIP'COOL

Hip compression cryotherapy belt

DESIGNED TO

- Relieving pain
- Reducing edema
- Limiting bleeding and hematoma
- Compatible with one day surgery
- Easy to use at home



Excell' ICE



Excell'ICE cold kits combine cryotherapy with modulated static compression using an inflation pump.

Easy to apply, its technology consists of alternating compressive cold during recovery and/or healing sessions, commonly recommended in the days following an injury or post-operation.

It is well known that external static compression with ice increases the extent of cooling at the surface of the skin.

A unique new concept where the compression chamber is directly integrated into the cold pack, making it much easier to use. This new technology makes it easier to carry out 20-30 minute cryotherapy-compressive recovery sessions. These portable devices are very practical to use at home or when travelling.

For athletes, this will ease musculoskeletal pain immediately after their matches or competitions and enable them to recover optimally, thereby reducing the time they spend out of competition due to pain, injury, fatigue, etc.



Recommended use:

20 min. by session every 4 hours
First session in Operating Room after dressing

Made in Europe

MD Medical Device

CE Compliant with European Union requirements

INDICATIONS

Post-operative, use in hip surgery.
Analgesic, anti-oedematous and anti-haemorrhagic effects.

CONTRAINDICATIONS

Allergy to cold, insensitivity to cold, Raynaud's syndrome or other vasospastic diseases, arteritis, cryoglobulinaemia, sickle cell disease, arteriosclerosis, diabetes, polyarthritis, skin sensitivity or hypercoagulation disorders, peripheral vascular diseases causing ischaemia or poor local blood circulation.

Belt

The universal hip belt may be required and indicated by surgeons for recovery after surgery or trauma. Indeed, with its adjustable immobilization abdominal and thigh straps, it allows the cold pack to be positioned according to the approach route. Effortless to put on, it offers a comfortable fit for the patient.

Hip compression cryotherapy device

Introducing our hip compression cryotherapy a cooling solution that is incredibly easy to apply inside the belt, the cold pack is compartmentalised to distribute the cold evenly and remains flexible to adapt to all body shapes. This device includes one single cold pack surrounding the hip, providing comprehensive coverage. Combine the cold therapy with compression in a single, convenient piece and achieve optimal compression using the included manual pump.

Excell' ICE



KHip'Cool (Abdominal / hip support belt and Compressive Cryotherapy Pack)

Ref. **ORT40HAN-ABD**



abdominal support belt



belt extension



thigh strap



thigh strap extension



Cold pack



Manual pump



Freezer bag



Suitcase

 **Orthonov**

3 rue de Villebon - 91160 Saulx-les-Chartreux - France
Phone: +33 1 60 13 65 73 - Email: contact@orthonov.com
SIRET: 820 720 068 00023 / VAT: FR01820720068

 Distributed by :