

# Universal Elbow Compressive Cryotherapy Device

## DESIGNED TO

- Immobilizing the elbow
- Fighting against pain
- Limiting bleeding and hematoma
- Managing post operative edema
- Easy to use at home



Excell'ICE



Excell'ICE cold kits combine cryotherapy with modulated static compression using an inflation pump.

Easy to apply, its technology consists of alternating compressive cold during recovery and/or healing sessions, commonly recommended in the days following an injury or post-operation.

It is well known that external static compression with ice increases the extent of cooling at the surface of the skin.

A unique new concept where the compression chamber is directly integrated into the cold pack, making it much easier to use. This new technology makes it easier to carry out 20-30 minute cryotherapy-compressive recovery sessions. These portable devices are very practical to use at home or when travelling.

For athletes, this will ease musculoskeletal pain immediately after their matches and competitions and enable them to recover optimally, thereby reducing the time they spend out of competition due to pain, injury, fatigue, etc.



### Recommended use:

20 min. by session every 4 hours  
First session in Operating Room after dressing

Made in Europe

**MD** Medical  
Device

**CE** Compliant with European  
Union requirements

## INDICATIONS

Management of lateral or medial epicondylalgia. After rehabilitation session. Analgesic, anti-oedematous and anti-haemorrhagic effects.

## CONTRAINDICATIONS

- Allergy to cold
- Insensitivity to cold
- Raynaud's syndrome
- Cryoglobulinemia
- Arteritis
- Diabetes
- Polyarthritis...

*Excell' ICE*

## Today's most user-friendly system: "Cold Pack" and "Compression Chamber" combined in a single device.

- Compression multiplies the effects of cold
- Cold packs remain flexible to adapt to body shape even at -4°F (-20°C)
- No frost to avoid skin burns
- Hypoallergenic Lycra with tightened mesh (included)
- 360° coverage around the Elbow in a single cold pack
- Freezer bag with zip slider for easy opening/closing (included)

## Elbow Cryo kit

Introducing the elbow Cryo Kit, a cooling solution designed for effortless application. Featuring a universal one-size-fits-all design, it conforms perfectly to any elbow. This kit is built with one single cold pack surrounding 360° around the elbow ensuring complete coverage. Its innovative design integrates cold therapy with compression, offering dual benefits in one convenient package. Utilize the included manual pump to achieve optimal compression and accelerate recovery.



## Universal Elbow Compressive Cryotherapy Device

Ref. ORT40PACK-C



Cold pack



Manual pump



Lycra



Freezer bag



Suitcase