

Wrist Brace + Compressive Cryotherapy Device

CRYOTHERAPY BENEFITS

- Limiting bleeding and hematoma
- Relieving pain
- Reducing edema
- Compatible with one day surgery
- Easy to use at home



Excell'ICE



Excell'ICE cold kits combine cryotherapy with modulated static compression using an inflation pump.

Easy to apply, its technology consists of applying compressive cold during recovery and/or healing sessions, commonly recommended in the days following an injury or post-operation.

It is well known that external static compression with ice increases the effectiveness of the cold to be diffused deeper into the joint.

A unique new concept where the compression chamber is directly integrated into the cold pack, making it much easier to use. This new technology makes it easier to carry out 20-30 minute cryotherapy-compressive recovery sessions. These portable devices are very practical to use at home or when travelling.

For athletes, this will ease musculoskeletal pain immediately after their matches or competitions and enable them to recover optimally, thereby reducing the time they spend out of competition due to pain, injury, fatigue, etc.



Recommended use:

20 min. by session every 4 hours
First session in Operating Room after dressing

Made in France

INDICATIONS

Post-operative care for wrist surgery, non-operated traumatology (scaphoid fracture, fracture of the base of one or more metacarpals of the long fingers, wrist sprain), after rehabilitation sessions.

Analgesic, anti-oedematous and anti-haemorrhagic effects.

CONTRAINDICATIONS

Cold allergy, cold insensitivity, Raynaud's syndrome, arteritis, cryoglobulinaemia, diabetes, polyarthritis, skin sensitivity disorders...

Excell' ICE

Today's most user-friendly system: "Cold Pack" and "Compression Chamber" combined in a single device.

- Compression multiplies the effects of cold
- Cold packs remain flexible to adapt to body shape even at -4°F (-20°C)
- No frost to avoid skin burns
- Hypoallergenic Lycra with tightened mesh (included)
- 360° coverage around the Wrist in a single cold pack
- Freezer bag with zip slider for easy opening/closing (included)



1 Wrist Cryo Kit

Introducing the Wrist Cryo Kit, a wrist cooling solution designed for effortless application. Featuring a universal one-size-fits-all design, it conforms perfectly to any wrist. This kit is built with one single cold pack surrounding 360° around the wrist ensuring complete coverage for your wrist. Its innovative design integrates cold therapy with compression, offering dual benefits in one convenient package. Utilize the included manual pump to achieve optimal compression and accelerate recovery.

2 Wrist Brace

This universal one-size wrist brace fits all and allows easy put on and provides great stability of the wrist. Thank to its adjustable straps it provides the necessary stability of the Wrist. Recommended for use post-wrist surgery or following injury, it aids in stabilization and promotes proper healing.

Wrist Brace + Compressive Cryotherapy Device

Ref. **ORT20POI-U**



Cold pack



Wrist Brace



Manual pump



Freezer bag



Lycra